


## FREQUENTLY ASKED QUESTIONS

# AMP MD SYSTEM

## USAGE & MAINTENANCE INSTRUCTIONS

---

Use roller after cleansing and toning your face in the evening. Use on face and/or neck at least three times a week. Gently roll the skin with 4-10 passes, changing direction with each pass. Total rolling time should be no more than one minute. For best results, use the roller on one section of your face at a time (forehead, left side, right side, chin and neck), completing each one before moving on to the next.



After each use, rinse the AMP MD Derma-Roller with very warm running water. Twice per week, clean thoroughly using AMP MD Purification Tablets, using a new tablet for each use:

- Fill Cleansing Container with very warm (not hot) water to marked fill line and drop one Purification Tablet into water.
- Securely lock the AMP MD Derma-Roller handle into the cap of the Cleansing Container and immediately immerse the AMP MD Derma-Roller into the cleansing solution.



- When the blue color fades, after approximately 15 minutes, remove AMP MD Derma-Roller and rinse thoroughly under very warm running water. Do not allow AMP MD Derma-Roller to remain in cleansing solution overnight.

**Product Storage:** Always store AMP MD Derma-Roller with the protective cap in place or locked into the dual-purpose Cleansing and Storage Container. Do not place the AMP MD Derma-Roller directly on hard surfaces or counter tops without the protective cap.

**For Your Safety:** AMP MD Derma-Roller should not be used on the eyes, on irritated or infected skin, and is not recommended for use on sensitive skin. Keep out of reach of children.